

Maximum Versatility
with Minimal Access



Mini-Thoracoscopy Set
for single-port medical thoracoscopy

Mini-Thoracoscopy Set

for single-port medical thoracoscopy

Richard Wolf Medical Instruments' unique Mini-Thoracoscopy Set enables minimally invasive diagnostics and treatment through one small incision—bringing the benefits of single-port, video-assisted interventions to the pulmonary space.

Benefits of the Mini-Thoracoscopy Set:

- Single piece, 0° operative endoscope with a 3.5 mm built-in working channel for single-port access and interventions
- Slim 5.5 mm outer diameter for a smaller incision
- Comprehensive set of compatible instrumentation for expanded treatment options
- Ergonomic endoscope and hand instrument design for enhanced comfort and control



| Description | Product No. |
|--|-------------|
| Endoscope | |
| Operative Telescope 5.5 mm w/ 3.5 mm Channel Bndl Includes: Op-endoscope 0°, Ø 5.5 mm, SL 215 mm (8920.401), automatic valve, i.d. 5.7 mm (8920.311), sealing membrane, Ø 17 mm (89.103), sealing cap, cap. 2.4-3.4 mm (89.01) | 89204015 |
| Trocar Sleeves and Trocars | |
| Trocar Sleeve (Flexible), straight, cap. 5.7 mm, WL 60 mm | 8906.051 |
| Trocar (Threaded), tip dull, sz. 5.5 mm, WL 83 mm | 8906.151 |
| Trocar Sleeve, straight, cap. 5.5 mm, WL 62 mm | 8919.333 |
| Trocar, tip dull, sz. 5.5 mm, WL 104 mm | 8919.3311 |
| Instruments and Accessories | |
| Double Spoon Forceps, Mono, Ø 3.5 mm Bndl Includes: Double spoon forceps insert, Ø 3.5 mm (8391216), sheath tube, insulated, Ø 3.5 mm, SL 330 mm (8391933), handle, pistol shaped, mono (83930074) | 83912167 |
| Dissection Forceps, Mono, Ø 3.5 mm Bndl Includes: Dissection forceps insert, Ø 3.5 mm (8391222), sheath tube, insulated, Ø 3.5 mm, SL 330 mm (8391933), handle, pistol shaped, mono (83930074) | 83912227 |
| Hook Electrode, mono, Ø 3.5 mm, WL 310 mm | 8379.452 |
| Coag Electrode (Button), mono, Ø 3.5 mm, WL 310 mm | 8379.462 |
| Suction Tube, Ø 2.8 mm, WL 450 mm | 8380.68 |
| Probe, Ø 3.5 mm, WL 250 mm | 8379.672 |
| Fiber Light Cable Bndl (Ø 2.5 mm) Includes: Fiber light cable, Ø 2.5 mm, TL 3 M (80662530), adapter projector side (8095.07), adapter endoscope side (809509) | 806625301 |

GUARDIAN
Protecting Performance



SERVICE BY RICHARD WOLF