

## Single piece instruments for increased procedural versatility



**5 Fr. Single Piece  
Semi-Rigid Instruments**

## 5 Fr. Single Piece Semi-Rigid Instruments

### Richard Wolf 5 Fr. Single Piece Semi-Rigid Instruments feature:

- A single piece design that requires no assembly prior to use
- An integrated flush port that allows for more thorough reprocessing of the instrument
- Six unique jaw patterns that accommodate a wide variety of hysteroscopic procedural applications
- Validation for steam sterilization
- The ability for disinfection using high-level disinfectant (HLD)



Jaw Pattern	Description	Outer Diameter Ø	Working Length (WL)	Product No.
	Grasping Forceps, Double-Action (Alligator)	5 Fr.	340 mm	<b>86422436</b>
	Biopsy Forceps, Double-Action (Cup)	5 Fr.	340 mm	<b>86422430</b>
	Tenaculum, Double-Action	5 Fr.	340 mm	<b>86422341</b>
	Scissors, Single-Action (Pointed)	5 Fr.	340 mm	<b>86422335</b>
	Scissors, Single-Action (Blunt)	5 Fr.	340 mm	<b>86422434</b>
	Scissors, Single-Action (Hook)	5 Fr.	340 mm	<b>86422416</b>