

All-in-one.




LED 50 HD II
Office Camera System

LED 50 HD II Office Camera System

All-in-one.

The agile LED 50 HD II Office Camera System is Richard Wolf's visualization solution for the office, featuring an all-in-one unit comprising the camera controller, light source and integrated recording system which complements an expansive portfolio of existing accessory components and endoscopic instrumentation.

With the LED 50 HD II Office Camera System, you can expect a space-saving, economical system with brilliant, high quality endoscopic imaging.



**21.5" HD
Medical
Monitor**

LED 50 HD II Camera Head

LED 50 HD II All-In-One Unit

UP-27MD Printer

**Condor Multi-Use Equipment Cart
with Condor Cart Video Kit**

LED 50 HD II Office Camera System

All-in-one.



spirit of excellence



LED 50 HD II All-In-One Unit

- All-in-one-unit combining camera controller, light source and integrated recording system allowing for a more economical and space-saving setup optimized for the office
- Plug-and-play for simple setup and intuitive operation
- High quality HD image (1920 x 1080 pixels)
- Bright and efficient LED lamp for homogenous, harmonized illumination
- Integrated recording system for storage of still images and videos on a USB stick

LED 50 HD II Camera Head

- 1/2" CMOS sensor with HD resolution (1920 x 1080 pixels)
- Ergonomic, lightweight design (0.39 lbs)
- Three programmable camera head buttons for quick and easy access to three functions
- Universal C-Mount threading for access to multitude of coupler options



LED 50 HD II Office Camera System

Ordering Information



spirit of excellence

All-In-One Unit

LED 50 HD II.....862-2519
Comprising: Camera Controller, LED Light Source,
Integrated Recording System

Camera Head

LED 50 HD II Camera Head.....740-2021

Couplers

RIWO Fixed Lenses, C-Mount

14 mm.....85261144

24 mm.....85261244

RIWO Angular Lenses, C-Mount

14 mm.....5257141

17 mm.....5257171

RIWO Zoom Lens, C-Mount

13-29 mm.....5261504

ENDOzoom Coupler, C-Mount

16-34 mm.....87-10-2020

Printer

Medical Grade Printer.....UP-27MD

Color Printing Pack.....5365.0211

Monitor

21.5" HD Medical Monitor.....FS-E2101D

Monitor Stand

Monitor Stand.....SU-600MD

Cart

Condor Multi-Use Equipment Cart.....2208030

Condor Cart Video Kit*.....2208029

Comprising: x1 Shelf, x1 VESA Monitor Mount, x1 Hardware Kit

*Please be advised that in order for the Condor Multi-Use Equipment Cart (PN: 2208030) to be compatible with a monitor, the Condor Cart Video Kit is needed (PN: 2208029).

