


Efficiency in a new size.



**Ultra-Mini PCNL
(UMP) Set**

Ultra-Mini PCNL (UMP) Set

Efficiency in a new size.



spirit of excellence

Ideal for effective and patient-friendly kidney stone treatment.

- High stone-free rate with low pressure and integrated, high-flow irrigation sheath.
- Small, 4 mm incision.
- Gentle and atraumatic intervention—the small puncture causes minimal dilation of the tissue.
- Lightweight design.
- Specially designed, state-of-the-art access sheath with water jet channel for creating turbulence of fragments and high-pressure.
- Excellent image quality and outstanding color reproduction.



Ordering Information

Product Number	Material Description
41.0665a	UMP Bundle Bundle includes: Ultra-Mini Nephroscope, all sheaths and obturators, cleaning basket, cleaning brush set, and sealing caps
41.0666a	UMP–Optic Irrigation Tube
41.0667a	UMP–Ultra-Mini Nephroscope 3.0 Fr., autoclavable, 0° direction of view, 85° field of view
41.0668a	UMP–6 Fr., Inner Sheath 220 mm working length
41.0670a	UMP–Water Jet Sheath, 11 Fr. 150 mm working length
41.0671a	UMP–Obturator for Water Jet Sheath, 11 Fr. 150 mm working length
41.0669a	UMP–7.5 Fr., Inner Sheath 220 mm working length
41.0672a	UMP–Water Jet Sheath, Standard, 13 Fr. 150 mm working length
41.0673a	UMP–Obturator for Water Jet Sheath, Standard, 13 Fr. 150 mm working length
41.0674a	UMP–Water Jet Sheath, Long, 13 Fr. 180 mm working length
41.0675a	UMP–Obturator for Water Jet Sheath, Long, 13 Fr. 180 mm working length
41.0676a	UMP–Cleaning Basket
41.0677a	UMP–Cleaning Brush Set
41.0678a	UMP–Sealing Cap for Laser Channel

To arrange a trial, contact your local Richard Wolf representative or call Customer Service at (800) 323-9653.