

All-in-one.



LED 50 HD II

Camera control unit

LED 50 HD II

Camera control unit

The agile LED 50 HD II camera control unit is Richard Wolf's visualization solution for the office, featuring an all-in-one unit comprised of a camera controller, light source and integrated recording system which complements an expansive portfolio of existing accessory components and endoscopic instrumentation.

With the LED 50 HD II camera control unit, you can expect a space-saving, economical platform with brilliant, high quality endoscopic imaging.



LED 50 HD II camera control unit

- All-in-one unit combining camera controller, light source and integrated recording system allowing for a more economical and space-saving setup optimized for the office
- Plug-and-play for simple setup and intuitive operation
- High quality HD image (1920 x 1080 pixels)
- Bright and efficient LED lamp for homogenous, harmonized illumination
- Integrated recording system for storage of still images and videos on a USB stick

Ordering Information:

LED 50 HD II

Comprising: Camera Controller, LED 50 HD II Camera Head, power cable, HDMI - DVI cable, USB flash drive 32 GB **862-2512**

Couplers

RIWO Fixed Lenses, C-Mount

14 mm **85261144**
24 mm **85261244**

RIWO Angular Lenses, C-Mount

14 mm **5257141**
17 mm **5257171**

RIWO Zoom Lens, C-Mount

13-29 mm **5261504**

LED 50 HD II camera head

- 1/2" CMOS sensor with HD resolution (1920 x 1080 pixels)
- Ergonomic, lightweight design (approx. 180 g)
- Three programmable camera head buttons for quick and easy access to three functions
- Universal C-Mount coupler for access to multitude of coupler options



Light Cables

2.5 mm Fiber Light Cable Bundle

Comprising: Adapter for Richard Wolf Light Source (8095.07), Adapter for Richard Wolf Endoscope (809509) and 2.5 mm Fiber Light Cable (80662530) **806625301**

Monitor

21.5" HD Medical Monitor **FS-E2102D**

